

# #690N/#691N Elite Care Cliners



ISO 9001-2000 CERTIFIED

## *Design, Comfort, Durability, and Value That Exceed Your Expectations*

Attention to detail and concern for comfort and quality – standards that Winco recliners are known for – place Winco Elite Care Cliners above the competition and make them among the best-selling chairs on the market today. Your Care Cliner now comes with LiquiCell®, an ultra-thin liquid-filled interface that aids in the process of preventing pressure sores without interfering with patient positioning or affecting stability.



#690N  
Elite Care Cliner



#691N  
XL Elite  
Care Cliner

For Use in Dialysis,  
Oncology, Acute Care,  
& Other Clinical Environments




Proudly Made in the USA!

# #690N/#691N Care Cliner




ISO 9001-2000 CERTIFIED

Adjustable elastic bands on back frame provide exceptional comfort.



Built-in, Ergonomic Headrest




Webbed Elastic and an Open Seat Base Provide Superior Comfort and Support.

Standard Trendelenberg Release on Back With Easy-Access Handle



Liquicell® Improves Blood Perfusion and Aids in Preventing Pressure Sores

Non-Porous, Fold-Down Side Tray Contains Small Spills  
Standard on Right Side  
Optional on Left Side



High-Impact, Nylon, 5" (12.7cm) Casters - Directional-Lock in Left Rear. All others are Total-Lock.

## VINYL SELECTION

(Printed vinyl colors may vary slightly from actual colors.)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Burgundy	Mauve	Blueridge	Grotto	Eucalyptus
<b>02</b>	<b>03</b>	<b>04</b>	<b>06</b>	<b>07</b>	<b>14</b>	<b>16</b>	<b>17</b>	<b>32</b>	<b>33</b>
<b>34</b>	<b>37</b>	<b>41</b>	<b>42</b>	<b>43</b>	<b>44</b>	<b>60</b>	<b>61</b>	<b>62</b>	<b>63</b>
Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress

## #690N/#691N SPECIFICATIONS

## OPTIONS

SEAT HEIGHT:	21" (53.3cm)
SEAT WIDTH:	21" (53.3cm) - #690N 25" (63.5cm) - #691N
SEAT DEPTH:	20" (50.8cm)
OVERALL HEIGHT:	48.5" (123.2cm)
OVERALL WIDTH:	34" (86.4cm)* - #690N 38" (96.5cm)* - #691N
BACK HEIGHT:	32.5" (82.6cm)
WEIGHT:	100 lbs. (45.4 kg) - #690N 123 lbs. (55.8 kg) - #691N
WEIGHT CAPACITY:	275 lbs. (124.7 kg) - #690N 450 lbs. (204kg) - #691N
FOOTPRINT:	34" x 41" (86.4cm x 104.1cm)* - #690N 38" x 41" (96.5cm x 104.1cm)* - #691N
FOOTPRINT (RECLINED):	34" x 76.5" (86.4cm x 194.3cm)* - #690N 38" x 76.5" (96.5cm x 194.3cm)* - #691N

\* ADD 10" (25.4cm) PER SIDE TABLE WHEN RAISED.




**TB-133**  
**Where Required**  
**#690NT #691NT**



**2nd Side Table on Left**  
**#690ND #691ND**



**Heat & Massage**  
**#690NH #691NH**  
**Heat ONLY**  
**#690NW #691NW**



**Headrest Cover**  
**#HC #HX**



**Square Leg Extensions**  
**#SE**



**IV Pole**  
**(factory installed)**  
**#690NIB #691NIB**