

#694N/#695N Dual Swing-Arm Elite Care Cliners



ISO 9001-2000 CERTIFIED

Design, Comfort, Durability, and Value That Exceed Your Expectations

Attention to detail and concern for comfort and quality - standards that Winco recliners are known for - place Winco Elite Care Cliners above the competition and among the best-selling chairs on the market today. Your Care Cliner now comes with LiquiCell®, an ultra-thin liquid-filled interface that aids in the process of preventing pressure sores without interfering with patient positioning or affecting stability.



#694N
Dual Swing-Arm
Elite Care Cliner



#695N
XL Dual
Swing-Arm
Elite
Care Cliner

For Use in Dialysis,
Oncology, Acute Care,
& Other Clinical Environments



Proudly Made in the USA!

#694N/#695N Dual Swing-Arm Elite Care Cliners




ISO 9001-2000 CERTIFIED

Standard Trendelenberg Release on Back With Easy-Access Handle




Built-in, Ergonomic Headrest

Webbed Elastic and an Open Seat Base Provide Superior Comfort and Support



Adjustable Elastic Bands on Back Frame Provide Exceptional Comfort



Liquicell® Improves Blood Perfusion and Aids in Preventing Pressure Sores

Non-Porous, Fold-Down Side Tray Contains Small Spills Standard on Right Side Optional on Left Side



Easy Swing-Arm Release



High-Impact, Nylon, 5" (12.7cm) Casters - Directional-Lock in Left Rear. All Others are Total-Lock.

VINYL SELECTION

(Printed vinyl colors may vary slightly from actual colors.)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Burgundy	Mauve	Blueridge	Grotto	Eucalyptus
02	03	04	06	07	14	16	17	32	33
34	37	41	42	43	44	60	61	62	63
Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress

#694N/#695N SPECIFICATIONS

SEAT HEIGHT:	21" (53.3cm)
SEAT WIDTH:	21" (53.3cm) - #694N 25" (63.5cm) - #695N
SEAT DEPTH:	20" (50.8cm)
OVERALL HEIGHT:	48.5" (123.2cm)
OVERALL WIDTH:	34" (86.4cm)* - #694N 38" (96.5cm)* - #695N
BACK HEIGHT:	32.5" (82.6cm)
WEIGHT:	104 lbs. (47.2 kg) - #694N 127 lbs. (57.6 kg) - #695N
WEIGHT CAPACITY:	275 lbs. (124.7 kg) - #694N 450 lbs. (204kg) - #695N
FOOTPRINT:	34" x 41" (86.4cm x 104.1cm)* - #694N 38" x 41" (96.5cm x 104.1cm)* - #695N
FOOTPRINT (RECLINED):	34" x 76.5" (86.4cm x 194.3cm)* - #694N 38" x 76.5" (96.5cm x 194.3cm)* - #695N

* ADD 10" (25.4cm) PER SIDE TABLE WHEN RAISED.

OPTIONS



TB-133
Where Required
#694NT #695NT



Heat & Massage
#694NH #695NH
Heat ONLY
#694NW #695NW



Square Leg Extensions
#SE



2nd Side Table On Left
#694ND #695ND



Headrest Cover
#HC #HX



IV Pole
(factory installed)
#694NIB #695NIB