

#840/#845 Task & Pediatric Tables



ISO 9001-2000 CERTIFIED



#840D
Pediatric Table
w/ Drawers



#840C



#840S



#840



#845C



#845S



#845

*Big Features
in a Little Table*



#845D
Task Table
w/ Drawers



Proudly Made in the USA!

#840/#845 Task & Pediatric Tables



ISO 9001-2000 CERTIFIED



VINYL SELECTION

(printed vinyl colors may vary slightly from actual color)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Black	White	Burgundy	Mauve	Blueridge	Grotto
02	03	04	06	07	08	09	14	16	17	32
OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
33	34	37	41	42	43	44	60	61	62	63
Eucalyptus	Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress

#840 SPECIFICATIONS

#845 SPECIFICATIONS

OPTIONS

OVERALL HEIGHT: 39" (99 cm)
 OVERALL WIDTH: 29.5" (75 cm)
 OVERALL LENGTH: 44" (112 cm)
 WEIGHT: Weight Varies Based on Options - Please Call
 WEIGHT CAPACITY: 400 LBS. (181.4 kg)
 FOOTPRINT: 29.5" x 44" (61 cm x 112 cm)

OVERALL HEIGHT: 36.5" (93 cm)
 OVERALL WIDTH: 29" (74 cm)
 OVERALL LENGTH: 42.5" (108 cm)
 WEIGHT: Weight Varies Based on Options - Please Call
 WEIGHT CAPACITY: 400 LBS. (181.4 kg)
 FOOTPRINT: 29" x 42.5" (74 cm x 108 cm)

Module Options

No Option Chosen - Plain Table

- S - Shelf (#840S or #845S)
- D - Drawers (#840D or #845D)
- C - Cabinets (#840C or #845C)