

#855/#855S Recovery Couches



ISO 9001-2000 CERTIFIED

Perfect Whenever and Wherever a Low-Height Table is Needed

Winco's recovery couches are a great choice for almost any medical facility or school nurse's office. With a height of only 19", but a weight capacity of 400lbs., practically any patient can sit or lie comfortably. The optional shelf keeps patient-accessible supplies nearby.



#855
Recovery Couch



#855S
Recovery Couch
w/ Shelf



Proudly Made in the USA!

#855/#855S Recovery Couches



ISO 9001-2000 CERTIFIED

Comfortable Medium Density Foam for Superior Comfort

Heavy Gauge Expanded Vinyl

High Strength Powder-Coated Steel Frame with a 400 lb. Weight Capacity



Adjustable Leveling Feet

Convenient Storage Shelf (#855S ONLY)

Standard Paper Holder (Paper Cutter Not Shown)

VINYL SELECTION

(printed vinyl colors may vary slightly from actual color)

Royal Blue	Taupe	Moss Green	Hunter Green	Grey	Black	White	Burgundy	Mauve	Blueridge	Grotto
02	03	04	06	07	08	09	14	16	17	32
OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL	OPTIONAL
33	34	37	41	42	43	44	60	61	62	63
Eucalyptus	Aubergine	Navy	Briar	Rosebud	Evergreen	Nightshade	Pyrite	Gypsum	Ice Blue	Cypress

#855/#855S SPECIFICATIONS

OVERALL HEIGHT: 19" (48cm)
 OVERALL WIDTH: 28" (71cm)
 OVERALL LENGTH: 72" (183cm)

WEIGHT: #855 - 73lbs. (33.1 kg)
 #855S - 98 lbs. (44.5 kg)

WEIGHT CAPACITY: 400 LBS. (181.4 kg)
 FOOTPRINT: 28" x 72" (71cm x 183cm)

OPTIONS



Shelf Option
855S