


Liberator And Liberator Hi Lo Illuminator Series Specifications



Letters reference dimensions shown in the right hand column.

Standard Features...

- Plexiglass viewing area and film retainer lift off for easy interior access.
- Film activated auto switches.
- Cluster switches and a master switch on 2 over 2 units and up.
- Gravity roller grip film retainer.
- Wolf's E-"Z" Bar™ for easy surface wall mounting. Simply mount the bar on the wall and slip the illuminator over it. See page 26 for details.

Light Intensities

Standard Liberator

- 4 bulbs per 14" x 17" viewing area.

Hi Lo Liberator

- Variable Light intensities (6 bulbs to 4 bulbs) per each 14" x 17" viewing area.

Optional Features

- Hospital grade plug.
- A frame for recessed mounting.
- A mobile stand.
- Custom room installations.

Ordering Information

- Use stock numbers to place your order
- If you want a recessed illuminator, order the appropriate Recessing Kit.

Description	Liberator Stock Numbers			Viewing Area	Wall Mounted Dimensions		
	Standard	Hi Lo	Recess Kit		A	B	C
Single Unit	22801	22901	21701	14" x 17" 35.6 x 43.2	22	14	3-1/4
2 Bank Unit	22802	22902	21702	28" x 17" 71.2 x 43.2	22	28	3-1/4
3 Bank Unit	22803	22903	21703	42" x 17" 106.8 x 43.2	22	42	3-1/4
4 Bank Unit	22804	22904	21704	56" x 17" 142.4 x 43.2	22	56	3-1/4
5 Bank Unit	22805	22905	21705	70" x 17" 178.0 x 43.2	22	70	3-1/4
6 Bank Unit	22806	22906	21706	84" x 17" 213.6 x 43.2	22	84	3-1/4
1 Over 1	22811	22911	21711	14" x 34" 35.6 x 86.4	44	14	3-1/4
2 Over 2	22822	22922	21722	28" x 34" 71.2 x 86.4	44	28	3-1/4
3 Over 3	22833	22933	21733	42" x 34" 106.8 x 86.4	44	42	3-1/4
4 Over 4	22844	22944	21744	56" x 34" 142.4 x 86.4	44	56	3-1/4
5 Over 5	22855	22955	21755	70" x 34" 178.0 x 86.4	44	70	3-1/4
6 Over 6	22866	22966	21766	84" x 34" 213.6 x 86.4	44	84	3-1/4



Wolf X-Ray Illuminators comply with all ACR Recommendations.